

SOFT-SOCK PROSTHETIC SOCK



- NEW!** Non-roll top!
- Looks great!
- Eliminates bunching!
- Patent pending.



SOFT-SOCK® PROSTHETIC SOCK

The original Soft-Sock... Now even better!

Available in X-STATIC® – the Silver Fiber® or original polyester fibers wick moisture away from skin and inhibit odor in the sock. Soft-Sock feels soft and cuddly worn next to the skin. New non-roll top looks great and eliminates bunching. Knit-in "Eye-Hole", reinforced distal toe option. Patent pending. Fewer sizes fit most regular sock wearers. (#A through #5 or #6 sock sizes). CARE TIPS: Wash Soft-Socks fleeced side in. No bleach. No liquid softeners. Tumble dry, cool temperatures.

CUDDLY 'N SOFT



3-D TOE!

Hemispherically knit three-dimensional distal shape for better fit.
US Patent No. 7,363,778

Knit-in "Eye-Hole", reinforced distal toe option. Patent pending.

LIGHT-WEIGHT	3-PLY	5-PLY	6-PLY	SIZE	FITS SOCK SIZE	LENGTH
SOFT-SOCK — POLYESTER/LYCRA®						
1SP1-CH-XS	—	—	—	Child X-Short	A	6" (15cm)
1SP1-CH-SH	1SP3-CH-SH	1SP5-CH-SH	1SP6-CH-SH	Child Short	A	8" (20cm)
1SP1-CH-12	1SP3-CH-12	1SP5-CH-12	—	Child Medium	A	12" (30cm)
1SP1-NA-SH	1SP3-NA-SH	1SP5-NA-SH	1SP6-NA-SH	Narrow Short	B, 0	10"–14" (25–35cm)
1SP1-NA-MD	1SP3-NA-MD	1SP5-NA-MD	1SP6-NA-MD	Narrow Medium	B, 0	16"–20" (40–50cm)
1SP1-NA-LG	1SP3-NA-LG	1SP5-NA-LG	1SP6-NA-LG	Narrow Long	B, 0	22"–28" (55–70cm)
1SP1-RG-XS	1SP3-RG-XS	1SP5-RG-XS	1SP6-RG-XS	Regular X-Short	1, 2	8"–12" (20–30cm)
1SP1-RG-SH	1SP3-RG-SH	1SP5-RG-SH	1SP6-RG-SH	Regular Short	1, 2	10"–14" (25–35cm)
1SP1-RG-MD	1SP3-RG-MD	1SP5-RG-MD	1SP6-RG-MD	Regular Medium	1, 2	16"–20" (40–50cm)
1SP1-RG-LG	1SP3-RG-LG	1SP5-RG-LG	1SP6-RG-LG	Regular Long	1, 2	22"–28" (55–70cm)
1SP1-WD-SH	1SP3-WD-SH	1SP5-WD-SH	1SP6-WD-SH	Wide Short	3 +	10"–14" (25–35cm)
1SP1-WD-MD	1SP3-WD-MD	1SP5-WD-MD	1SP6-WD-MD	Wide Medium	3 +	16"–20" (40–50cm)
1SP1-XW-SH	1SP3-XW-SH	1SP5-XW-SH	—	X-Wide Short	4, 5+	10"–14" (25–35cm)
1SP1-XW-MD	1SP3-XW-MD	1SP5-XW-MD	—	X-Wide Medium	4, 5+	16"–20" (40–50cm)

HOLE-IN-TOE SOFT-SOCK — POLYESTER/LYCRA®

1SPS-CH-XS	—	—	—	Child X-Short	A	6" (15cm)
1SPS-CH-SH	1SP3-SC-SH	1SP5-SC-SH	1SP6-SC-SH	Child Short	A	8" (20cm)
1SPS-NA-SH	1SP3-SN-SH	1SP5-SN-SH	1SP6-SN-SH	Narrow Short	B, 0	10"–14" (25–35cm)
1SPS-NA-MD	1SP3-SN-MD	1SP5-SN-MD	1SP6-SN-MD	Narrow Medium	B, 0	16"–20" (40–50cm)
1SPS-RG-XS	1SP3-SR-XS	1SP5-SR-XS	1SP6-SR-XS	Regular X-Short	1, 2	8"–12" (20–30cm)
1SPS-RG-SH	1SP3-SR-SH	1SP5-SR-SH	1SP6-SR-SH	Regular Short	1, 2	10"–14" (25–35cm)
1SPS-RG-MD	1SP3-SR-MD	1SP5-SR-MD	1SP6-SR-MD	Regular Medium	1, 2	16"–20" (40–50cm)
1SPS-RG-LG	1SP3-SR-LG	1SP5-SR-LG	1SP6-SR-LG	Regular Long	1, 2	22"–28" (55–70cm)
1SPS-WD-SH	1SP3-SW-SH	1SP5-SW-SH	1SP6-SW-SH	Wide Short	3 +	10"–14" (25–35cm)
1SPS-WD-MD	1SP3-SW-MD	1SP5-SW-MD	1SP6-SW-MD	Wide Medium	3 +	16"–20" (40–50cm)

SOFT-SOCK WITH X-STATIC® — POLYESTER/X-STATIC®/LYCRA®

1SX1-NA-SH	1SX3-NA-SH	1SX5-NA-SH	—	Narrow Short	B, 0	10"–14" (25–35cm)
1SX1-RG-XS	1SX3-RG-XS	1SX5-RG-XS	—	Regular X-Short	1, 2	8"–12" (20–30cm)
1SX1-RG-SH	1SX3-RG-SH	1SX5-RG-SH	—	Regular Short	1, 2	10"–14" (25–35cm)
1SX1-RG-MD	1SX3-RG-MD	1SX5-RG-MD	—	Regular Medium	1, 2	16"–20" (40–50cm)
1SX1-WD-MD	1SX3-WD-MD	1SX5-WD-MD	—	Wide Medium	3 +	16"–20" (40–50cm)

