

# POST-OP PROSTHETIC SOCKS AND SUPPLIES



**X-static** Made with 10% X-Static® Antimicrobial Technology.

## SILVER A/K AND B/K POST-OP SOCKS

Prosthetic interface for immediate post operative care.\*  
 Made with 10% X-Static® Antimicrobial Technology. 3rd party testing showed 99.98% reduction of bacteria directly on the surface of the sock (vs 93% for control sample)\*\*

X-Static® is made by permanently bonding 99.9% pure silver to the surface of fibers. Socks powered by X-Static® inhibit the growth of bacteria and fungi directly on the surface of the fabric. Because it is permanent, it will not wash out and performs for the life of the product.

Packed in individual polyethylene bags, four (4) to a box. Fleeced inside. Rolled for application. 3-ply weight. Latex free materials. 87.5% Polyester, 10% X-Static® Silver, 2.5% Lycra®

Patent Pending.

\*(Not a primary sterile dressing)

\*\* Data on file

	LENGTH	TOP WIDTH	TOE WIDTH
<b>1PL3-X4-18</b>	18" (45cm)	5" (12.5cm)	4" (10.5cm)
<b>1PL3-X5-18</b>	18" (45cm)	6" (15cm)	5" (12.5cm)
<b>1PL3-X6-18</b>	18" (45cm)	8" (20cm)	6" (15cm)
<b>1PL3-X7-18</b>	18" (45cm)	9" (22.5cm)	7" (17.5cm)



Post-op socks for immediate post operative prosthetic fittings

## ORIGINAL UNSTERILIZED A/K AND B/K POST-OP SOCKS

Packed in individual polyethylene bags, four (4) to a box. 3-Ply Weight. Latex free materials. 97% Acrylic, 3% Lycra. Modified box toe. Fleeced inside. Rolled for application.

	LENGTH	TOP WIDTH	TOE WIDTH
<b>1PL3-B5-14</b>	14" (35cm)	5" (12.5cm)	3.50" (9cm)
<b>1PL3-B5-18</b>	18" (45cm)	5" (12.5cm)	4" (10cm)
<b>1PL3-05-18</b>	18" (45cm)	6" (15cm)	5" (12.5cm)
<b>1PL3-26-18</b>	18" (45cm)	8" (20cm)	6" (15cm)
<b>1PL3-37-18</b>	18" (45cm)	9" (22.5cm)	7" (17.5cm)

NOTE: Socks knitted from this material are sized to fit per sizing of other materials such as wool. Due to the stretch and contraction of the material the flat measurements will be approximately 10% less in both length and width.

## SYMES AND CHOPART UNSTERILIZED POST-OP SOCKS

For immediate post operative use. 3-Ply acrylic/lycra, seamless for optimum fitting qualities. Warning — Must be gas autoclaved to be sterile. Latex free materials. 97% acrylic, 3% Lycra. Fleeced inside. Sold each.

SYMES	LENGTH	TOP WIDTH	ANKLE WIDTH	BALL
<b>1SL3-E1-24</b>	24" (60cm)	7.50" (19cm)	4.50" (11cm)	5" (12.5cm)
<b>1CL3-E1-24</b>	24" (60cm)	7" (18cm)	4" (10cm)	6" (15cm)

